



## **BPL ENDO VIEW**

Reusable Bronchoscope



*Happier Living Everyday*



Follow us on



[www.bplmedicaltechnologies.com](http://www.bplmedicaltechnologies.com)



# BPL ENDO VIEW

## Reusable Flexible Bronchoscope



### Display

Horizontal and vertical rotation angle:  $\geq 120^\circ$   
HD full view screen,  
6 function buttons, USB port,  
HDMI port for big display.



### Quick Handlings

1 lock design  
Fast installation and disassembling  
Stay firm and efficient



### Angulation Control

Flexible angulation control  
Precise location



### Scope

True 1:1 Torque Ratio for Easy operation



### Camera System

Photo and Video recording function,  
Easy for clinical teaching and training



### Suction Button

Fast Connection; Smooth Suction



### Handle

Ergonomic streamline design  
Anti-slip and improve hand held comfort



### Working Chanel

Inner diameter:  
1.2mm/2.6mm/3.0mm available

### Insertion part

Outer diameter: 2.5mm/3.9mm/5.2mm/5.8mm  
Fracture resistant, smooth, flexible and durable



### Distal tip part

Dual deflections  
New CMOS Technology  
High brightness LED light source  
160 up 130 down



## Reference

Size		M58	M52	M39	M28
Scope Body	Illumination	LED			
	Joint Outer Diameter	5.8mm	5.2mm	3.9mm	2.8mm
	Instrument Chanel Inner Diameter	3.0mm	2.6mm	1.2mm	0mm
	Working Lenth	650mm			
	Anti-Fog	Anti-fog when it's on, no need for preheating.			
Operating system 、 Software System upgrades are free for life					
Display	Size	3.5inch/8inch			
	Storage	8GB SD Card			
	Shoot	Support			
	Video Record	Support			
	White Balance	Support			
	Freeze/Backwards	Support			
	Video Synchronous Out put	Support			

**Note:** » M28/M39 is not applicable for endoscopic diagnosis and treatment accessories.

» M28/M39/M52/M58 is not suitable for high-frequency electrotherapy surgery and laser hemostasis.

» Continuous working time in more than two hours, but it depends on ambient temperature and battery condition.

# BPL ENDO VIEW

## Reusable Bronchoscope



### Advantages

#### Portable

Light weight, Easy to operate, Saving time for EMS, Convenient for transportation and storage

#### Multifunction

Available for the intubation of difficult airway, Inspection of airway and double chamber, Alveolar lavage, Biopsy, Treatment of pneumothorax, Stent Implantation, Removal of tracheal foreign bodies, Other treatments for the bronchi, Such as Local radiotherapy and chemotherapy, Medication Injection

#### Economic

One screen for four scopes, Four working channel sizes for option

#### Durable

Electronic camera system with LED, New technology for the photo processing: CMOS, Higher quality pictures, Last longer than traditional scope

#### Application

Applied to the departments such as: ICU, NICU, Respiration, Anesthesia, Emergency etc.

### Accessories available: Trolley



**ISO 13485:2016 CERTIFIED COMPANY**  
**ISO 9001:2015 CERTIFIED COMPANY**

**BPL MEDICAL TECHNOLOGIES PRIVATE LIMITED**  
Regd. Office: 11th KM, Bannerghatta Road,  
Arakere, Bangalore - 560076, India.  
Toll Free: 1800-4252355  
Website: [www.bplmedicaltechnologies.com](http://www.bplmedicaltechnologies.com)  
For Enquiries: [sales.medical@bplin](mailto:sales.medical@bplin)  
CIN: U33110KA2012PTC067282



*Happier Living Everyday*