

Medical Irrigation Pump



Technical Parameters

- Hysteroscopy intra-cavity pressure management
- Irrigation and suction pump inbuilt single console
- 2 Dedicated pressure sensor for irrigation and suction
- Flow management system with (100ml - 1000ml)
- Smart Micro Controller based
- disposable tubings D EHP-free.
- 4" touch screen with user defined controls
- input voltage - 85-240V~AC 50 Hz
- Dimensions:- WXHXD 155mm x 155mm x 310mm
- Weight: 6.5kg approx
- Flow settings - 0.1-1.0L/min
- Pressure settings - 5-400 mmHg
- Suction- 0.65 bar
- User defined programable control
- Arthroscopy - user defined
- Urology - user defined
- Hysteroscopy - user defined

MEDITRAX EQUIPMENTS proud to introduce the state of art MEDICAL IRRIGATION /SUCTION PUMP for Hystroscopy / Arthroscopy / Laparoscopy/ Urology. 2 years of our hard work and research from our strong R&D team, we have developed this product.

it contains user defined programable controls which gives the user provision to work on their best settings according to their case, easy to program modes just set the desired pressure and flow then hold the mode for 5 seconds to save the settings.

