



India Corporate Office

Beurer India Private Limited
523-524, JMD megapolis, Sector- 48, Sohna Road,
Gurugram, Haryana, India - 122018

Beurer Headquarters

Beurer GmbH
Söflinger Str. 218 · 89077 Ulm
Germany



TOLL FREE NO
1800-208-8284

www.beurerindia.com

Blood Pressure

What is blood pressure?

For your blood to reach all your body organs, it needs to be pumped from the heart at a certain level of pressure so that the blood circulation in your body is maintained.

- Systolic pressure-The pressure exerted when the heart contracts and pumps blood.
- Diastolic pressure - The pressure at the moment the heart relaxes.

What is the need to measure Blood Pressure?

- Every 5th adult over 40 years old is affected by high blood pressure.
- High blood pressure is a risk factor that can lead to the development of severe illnesses. The normalization of blood pressure prevents strokes, lowers the risk of cardiac arrest, and reduces the risks of memory loss.
- By means of regular measurement, elevated blood pressure can be detected in time and treated.

Basic rules for measurement

- 1- Measure under the same conditions and at approx. the same time of the day to allow comparability.
- 2- Do not take a measurement within 30 minutes after eating, drinking, smoking, or exercising. Perform the measurement in absolute peace (mental, physical) after 5 - 10 minutes of relaxing.
- 3- Sit comfortably or lay down, do not stand.
- 4- Hold arm still (e.g., rest it on the table)
- 5- Do not speak, clench fists, or move.
- 6- Just relax before and during measurement



Upper Arm Blood Pressure Monitors



MRP: ₹2240/-

BM 27

Upper Arm Blood Pressure Monitor

Beurer BM 27 upper arm blood pressure monitor makes it easier to perform blood pressure check at home with an advanced data tracking feature. In an ergonomic design, this blood pressure measuring device is laced with cool features,


- Easy-to-read display
- Saves and provides the average of all measurements.
- 4 x 30 memory spaces
- Arrhythmia detection and risk indicator
- Date and time / automatic switch-off

BM 28

Upper Arm Blood Pressure Monitor

BM 28 blood pressure monitor makes it convenient to measure blood pressure in a domestic environment. This device comes with both battery as well as mains operation and offers various other features,

- Average of all recorded measurement
- A readable display
- 4x30 memory space
- Arrhythmia detection and risk indicator




MRP: ₹2760/- (with adaptor)

MRP: ₹2380/- (without adaptor)



MRP: ₹2360/-



NEW LAUNCH

MRP: ₹2095/-

BM 29

Upper Arm Blood Pressure Monitor

One for all- BM 29 comes in multicolor.

The universal cuff means the device is also suitable for larger upper arm circumferences.

- Easy to read display.
- 4x30 memory space
- Arrhythmia detection and risk indicator
- Date and time / automatic switch-off

BM 46


Upper Arm Blood Pressure Monitor

BM 46 is a fully automatic blood pressure monitor having an XL display.

- Clear, easy to read display
- Gives average values of measurement
- 4x30 memory space
- Arrhythmia detection and risk indicator

Features	BM 27	BM 28	BM 29	BM 46
Deflation Technology	✓	✓	✓	✓
Memory space: 4x30	✓	✓	✓	✓
Device Colour- white	✓	✓	X	✓
Mains adaptor	✓	✓	✓	X
Arrhythmia & Risk indicator	✓	✓	✓	✓
Cuff size 22-36 cm	✓	✓	✓	✓
Resting Indicator (HSD)	X	✓	✓	X
Storage Box	X	X	X	✓
Batteries	4x 1.5 V AA	4x 1.5 V AA	4x 1.5 V AA	4x 1.5 V AAA





MRP: ₹5610/-

BM 54 *Bluetooth®*

With the BM 54 upper arm blood pressure monitor, measurement values can be transferred to the smartphone via Bluetooth using the beurer health manager pro app. Other features include-

- Universal cuff, including for larger upper arm circumferences up to 44 cm
- Risk indicator: Classification of measurements with coloured scale
- Wireless transfer of measurements
- Clear, easy-to-read XL display
- 2x60 memory space
- Arrhythmia detection



MRP: ₹11250/-

BM 81 *Easy Lock*

The innovative easy Lock system makes it easy for you to fix the device on your upper arm with the perfect fit and secure hold – without the need for a hook-and-loop fastener. The measured values can be conveniently transferred to the "beurer Health Manager Pro" app and monitored long-term.

- Innovative easy Lock cuff: put it on, tighten and start.
- No hoses, no cables – easy and precise measurement.
- Inflation technology LED risk indicator: classification of measurements with coloured LED scale.
- 2x120 memory space.

BM 96 *with ECG*

With the multi-functional BM 96 device, you can keep track of your heart health. In addition to the fully automatic blood pressure and pulse measurement on the upper arm, your heart rhythm is also recorded. The “beurer Health Manager” app is the optimal way to monitor your values.

- 1-channel ECG device for recording heart rhythm
- Advanced arrhythmia detection and risk indicator
- Resting indicator for precise measurement results
- Easy operation via 4 sensor touch buttons
- 2x60 memory space
- Easy-to-read XL display
- Inflation technology



MRP: ₹12850/-

Features	BM 54 <i>Bluetooth®</i>	BM 81 <i>Easy lock</i>	BM 96 <i>with ECG</i>
Deflation Technology	✓	<i>Inflation technology</i>	<i>Inflation technology</i>
Memory space:2x60	✓	2x120	✓
Software/app Health manager	✓	✓	✓
Batteries- 4x 1.5V AAA	✓	✓	✓
LED Risk Indicator	X	✓	✓
Arrhythmia	✓	✓	✓
Illum. display	X	✓	✓
Cuff size 22-42cm	✓	24 – 40 cm	✓

Wrist Blood Pressure Monitors



MRP: ₹1999/-



MRP: ₹2380/-

BC 28

Wrist Blood Pressure Monitor

Beurer BC 28 wrist blood pressure monitor is straightforward and easy to use. These handy devices prove themselves to be particularly reliable companions on the move, meaning you can monitor your blood pressure regularly. Fully automatic blood pressure and pulse measurement on the wrist.

- Average of all saved measurements
- Risk indicator & arrhythmia detection
- For wrist circumferences from 14.0 to 19.5 cm
- 2 x 60 memory spaces
- Date and time/automatic switch-off
- Low battery indicator

BC 30

Wrist Blood Pressure Monitor

This wrist blood pressure monitor is a fully automatic device that is developed to capture the blood pressure and pulse measurements on your wrist. It has an easy-to-read screen that displays the readings along with the date and time. The device also has a Velcro strap fastener for efficient use.

- Average of the last 3 measurements
- 2 x 60 memory spaces
- Risk indicator and arrhythmia detection
- For wrist circumferences from 13.5 to 19.5 cm
- Date and time/automatic switch-off
- Low battery indicator



MRP: ₹2380/-

BC 32

Wrist Blood Pressure Monitor

Beurer BC 32 Wrist Blood Pressure Monitor (Black) has been especially designed for keeping an easy check on your blood pressure right from the vicinity of your home. The digital screen display makes the measurement an easy task.

- Average of all saved measured values
- Average of morning and evening blood pressure over the last 7 days
- 2 x 60 memory spaces
- Risk indicator and arrhythmia detection
- For wrist circumferences from 13.5 to 19.5 cm

Features	BC 28	BC 30	BC 32
Deflation Technology	✓	✓	✓
Memory space: 2x60	✓	✓	✓
Cuff size 14-19.5cm	14.0 – 19.5 cm	✓	✓
Risk indicator & arrhythmia	✓	✓	✓
Illuminated display	X	X	X
Batteries 2x 1.5V AAA	✓	✓	✓
Storage box	✓	✓	✓
Displays all saved values & the last 7 mornings and evenings	✓	last 3 measurements	✓

Body Temperature

What is fever?

Fever is the body defence mechanism against infections. Its measurement is dependent on the measuring point standard temperature spans measured at several locations.

Forehead temperature: 35.8°C to 37.6°C

Ear temperature: 36.0°C to 37.8°C

Oral temperature: 36.0°C to 37.4°C

Rectal temperature: 36.3°C to 37.8°C

Temperature is affected by:


- 1- Metabolism/age (small children have higher & more variable temperatures, this regulates with age)
- 2- Clothing in combination with outside temperature
- 3- Time of day (lowest in the morning and increases throughout the day)
- 4- Physical activity
- 5- Measurement point



Thermometers



MRP: ₹255/-



MRP: ₹500/-

FT09/1

Clinical Thermometer

This digital thermometer offers high precision and accuracy. It does not contain mercury and is completely waterproof. After immediate usage it can be easily disinfected. The thermometer also comes with a protective cover which prevents it from harmful germs and bacteria. When not in use it gets automatic switched off.

- Contact-measurement technology
- Display in °F.
- Waterproof
- No mercury, without glass
- Can be disinfected.
- 1 memory space
- Automatic switch-off

FT 13

Flexi Thermometer

The clinical thermometer has an extra-flexible measuring tip. The temperature alarm with warning tone from 37.8°C is practical and makes assessing the measurements easier.

- Extra-flexible measuring tip
- Measuring time approx. 30 sec.
- Contact-measurement technology
- Temperature alarm: warning tone from 37.8°C
- Displays measurements in °C and °F
- Waterproof
- No mercury, without glass
- Can be disinfected
- 1 memory space
- Automatic switch-off
- Acoustic signal when measurement is complete.



MRP: ₹2400/-

FT 65

Multifunctional Thermometer

Beurer Multi-Functional Infrared Thermometer is a smart medical device that has been developed to measure the user's forehead temperature or the surface temperature of objects and liquids in seconds. Additionally, it comes with a large, easy-to-read digital screen to display the readings.

- 6-in-1 function: ear, forehead and surface temperature, temperature alarm, date, and time, 10 memory spaces
- Up to 37.9°C = green, from 38.0°C = red
- Measurement in seconds
- Infrared measurement
- LED temperature alarm from 38°C
- Displays measurements in °C and °F.
- No mercury, without glass
- 10 memory spaces with date and time



MRP: ₹2500/-

FT 85

Non-Contact Thermometer

The Beurer FT 85 is a contact less infrared clinical thermometer that can be used to measure the body temperature of any individual without direct skin contact. This thermometer can also be used to detect the temperature of objects and your surroundings with absolute precision and accuracy.

- Measurement of body, ambient and surface temperature
- Hygienic and safe – measurement without skin contact
- Measurement in seconds
- Infrared measurement
- Temperature alarm: from 38°C
- Displays measurements in °C and °F.
- 60 memory spaces (body temperature)

Features	FT09/1	FT 13	FT65	FT85
Measurement Range	35.5°C – 42.0°C	35.5°C – 42.0°C	35°C – 42°C	34°C – 42.2°C
Accuracy	± 0.2°C	± 0.2°C	± 0.2°C	±0.2°C
Date and time	X	X	✓	✓
Memory	1 (last value)	1 (last value)	10	60
Displays both °C & °F	only °F	✓	✓	✓



Respiratory System

Respiratory infections are prevalent all around the world. Unfortunately, none of us are immune to it, especially when flu season starts.

Two methods to tackle lung infections are **Steam Inhalation** and **Nebulizer**.

What Is Steam Inhalation?

Steam inhalation is a prevalent home remedy to achieve relief from cold and sinus infections. It involves inhaling steam to open the nasal passages. The warm air loosens the mucus in your nasal pathways and throat. By clearing the passages, it provides you relief from symptoms of swollen blood vessels.

What Is Nebulizer?

A nebulizer converts liquid medicines into a mist. All you have to do is wear the oxygen mask and breathe in the mist. It can treat various lung conditions such as chest tightness and wheezing.

Most doctors prescribe nebulizers for asthma patients. But that is not the only use of the machine. It also serves as breathing therapy for many. .

Particle Size

The general difference between the two technologies is that the tiny respirable particles of the active substances which are continuously produced during use, reach different parts of the airways.

- $>6\text{ }\mu\text{m}$: particle only stays in mouth, throat, and nose.
- $<6\text{ }\mu\text{m}$: particle goes to deeper pulmonary wings.

↳ **All Beurer nebulizers reach $< 6\text{ }\mu\text{m}$**



Oxygen Saturation is a measure of how much oxygen is in your blood. It is a crucial measure of how well the lungs are working.

For most healthy adults, a normal oxygen saturation level is between 95% and 100%. An "O2 sat" level below this range requires medical attention because it means the body isn't getting enough oxygen.

Respiratory System

Oxygen Saturation

Pulse Oximeters


The pulse oximeter, or Pulse Ox, is an electronic device that measures the saturation of oxygen carried in your red blood cells.

It can rapidly detect even small changes in oxygen levels. These levels show how efficiently blood is carrying oxygen to the extremities furthest from your heart, including your arms and legs.

The purpose of pulse oximeter is to see if your blood is well oxygenated.



Nebulizer



MRP: ₹1990/-



MRP: ₹3425/-

IH 18

Nebulizer

IH 18 is a compressor nebulizer that uses compressed air technology to nebulise liquid medication. It helps in the speedy and effective treatment of upper and lower respiratory tracts.

- Compressed air technology
- For treating the upper and lower airways, colds, asthma, respiratory diseases
- Nebulisation performance: **approx. 0.4 ml/min**
- Short inhalation time
- High proportion of respirable particles
- Particle size $3.16\mu\text{m}$
- Atomizer size 8 ml

IH 21

Nebulizer

Beurer IH21 is used to treat the upper and lower respiratory tract, colds and asthma, and respiratory diseases. This Beurer Nebulizer Machine offers a high rate of medication vaporization (approx. 0.3 ml/min), making the inhalation time shorter.

- Compressed air technology
- For treating the upper and lower airways, colds, asthma, respiratory diseases
- Nebulization performance: approx. 0.3 ml/min
- Working pressure: approx. 0.8 to 1.45 bar
- Extensive accessories
- Mains operation
- Particle size $3.07\mu\text{m}$



MRP: ₹5330/-

IH 58 Nebuliser

IH 58 Nebulizer is an essential medical device that is developed especially for individuals who suffer from cold and other respiratory illnesses. It has a high nebulization capacity of about 0.25ml/min. The device operates at a pressure ranging between 0.25-0.5 bar.

- Compressed air technology
- For treating the upper and lower airways, colds, asthma, respiratory diseases
- With mains operation
- Ideal when travelling
- Nebulisation performance: approx. 0.25 ml/min
- Working pressure: approx. 0.25 to 0.5 bar
- Incl. mains adapter with micro-USB cable

Features	IH 18	IH 21	IH 58
Compressed air Technology	✓	✓	✓
Mains operation	✓	✓	✓
Nebulisation performance	0.4 ml/min	0.3 ml/min	0.25 ml/min
Particle size	3.16 µm	3.07µm	4.12 µm
Operating pressure: 0.8-1.45 bar	✓	✓	0.25-0.5 bar
Baby/ children version	X	X	✓

FS 50

Facial sauna

FS50 Facial Sauna and Steam Inhaler naturally cleans and moisturizes your skin. Along with facial care, beurer sauna is also used to give relief from cold as it acts as a steamer for cold perfectly with its steam attachment. It helps in treating sinus or other cold effects. To boost up the effect, beurer FS50 also allows adding essential oils or scent in its scent container.

- Ideal for both cosmetic facial care and for inhalation with the appropriate steam attachment for mouth and nose inhalation.
- Opens pores for more effective cleansing.
- Also, work as a steamer for cold conditions.
- You can use it as a steamer for nose.
- Scent container for aroma application
- 2 output levels
- Mains operation



MRP: ₹3710/-

SI 40

Steam Vaporizer

This simple and easy to use steam vaporizer is intelligently designed for warm steam treatment during the colder seasons. The vaporizer is ultra-quiet and nebulizes liquids such as saline solution and essential oils.

- Soothing treatment during the cold season
- Includes flexible universal mask for the mouth and nose.
- Steam setting can be varied for individual treatments.
- Effective application due to constant temperature
- Short application time
- Simple one-button operation
- Mains operation



MRP: ₹4045/-

PO 30

Pulse Oximeter

PO 30 Pulse Oximeter is a lightweight device that is used to monitor the SpO2 level and pulse rate of an individual. This non-invasive tool attaches painlessly to your fingertip in order to offer accurate and precise readings. It is very convenient to use and comfortable to carry securely while travelling.


- For checking and regularly monitoring oxygen saturation
- Also suitable at higher altitudes
- Graphic pulse display
- Small and light
- Low battery indicator
- Incl. strap and belt pouch
- Colour display with 4 available views

PO 40


Pulse Oximeter

With this small pulse oximeter, one can easily determine the arterial oxygen saturation, heart rate and Perfusion Index. The device is small and light- for use at home and on the move.

- Measurement of SpO2 as a potential vital sign in acute illnesses, heart rate (pulse) and perfusion index (PI)
- Colour display with 7 available views
- Adjustable display brightness
- Symbol display for inaccurate measurement
- Particularly suitable for persons with: Heart failure, COPD, Bronchial asthma
- Low battery indicator
- Incl. strap, belt pouch



MRP: ₹2380/-



MRP: ₹3875/-

Features	PO 30	PO 40
Measurements Range: SpO2- 0%-100%	✓	70%-100%
Measurement Range: Pulse rate- 30-250 bpm	✓	✓
Displays Perfusion Index (PI)	X	✓
Available views- 4	✓	✓
Batteries- 2x1.5V AAA	✓	✓

Hearing Amplifiers

V/s

Hearing Aids

Hearing Amplifiers

Hearing amplifiers are wearable devices that use a microphone to pick up nearby sounds. They then process this sound and deliver it to a person's ear at a louder volume.

Hearing Aid

Hearing aids are devices that make sounds louder so that people with mild-to-severe hearing loss can hear and more fully participate in everyday activities.

A hearing aid has a microphone, amplifier, and speaker. The microphone picks up the sound, and the sound becomes electrical signals that get sent to an amplifier. The amplifier then sends the sound waves to the ear via the speaker.



Hearing Amplifier



MRP: ₹1735/-



MRP: ₹2500/-

HA 20

Hearing Amplifier

This hearing machine device is ideal for individuals with restricted hearing ability. It fits comfortably behind the ear. The hearing amplifier helps you hear louder and clearer again. The discrete earpiece amplifies all sounds both indoors and outdoors.

- Individual adjustment to the ear canal
- Ideal for restricted hearing ability
- Frequency range: 200 to 5000 Hz
- Amplification: max. 40 dB
- Maximum volume of 128 dB
- 3 attachments to individually adjust to the ear canal.

HA 50

Hearing Amplifier

This ergonomically designed hearing amplifier is best suited for people enduring restricted hearing ability. It comfortably fits behind your ear and is hardly visible. The low noise rendering technology amplifies all sounds indoors and outdoors for sharper hearing.

- Extended frequency range
- Individual adjustment to the ear canal
- Ideal for restricted hearing ability
- Ergonomic fit behind the ear, barely visible
- Low-noise rendering
- Frequency range: 100 to 6000 Hz
- Amplification: max. 40 dB
- Volume: max 128 dB



NEW LAUNCH

MRP: ₹7900/- Single

MRP: ₹16900/- Pair

HA 81

Re-chargeable Hearing Amplifier

The digital hearing amplifier set amplifies the volume of all sounds indoors and outdoors and is ideal if you are hard of hearing.

- Rechargeable hearing amplifiers with exchangeable batteries.
- Handy set of 2 for amplification on both sides
- 2 different hearing programmes
- With micro USB Charger
- 4 attachments each to individually adjust to the ear canal
- Frequency range: 200- 6700 Hz
- Amplification Max: 48 dB
- Volume: max 121 dB
- Charging time: 2 Hours

Features	HA 20	HA 50	HA 81
Type	RIC	RIC	BTE
Signal	Analog	Analog	Digital
Amplification	Max 40 dB	Max 40 dB	Max 48 dB
Max volume	128 dB	128 dB	121 dB
Frequency range	200- 5000 Hz	100-6000 Hz	200-6700 Hz
Batteries	2x 1.4V PR48	2x 1.4V PR48	2x / 4x Li-Ion re-chargeable batteries
Battery Life	75 hours	60 hours	≥ 26 hours
Attachments	3	3	4

Electro stimulation (TENS)

What is Electro stimulation?

Electro stimulation is the imitation of impulses in our bodies which are transferred through the skin to nerve and muscle fibres via electrodes.

TENS: Transcutaneous Electrical Nerve Stimulation => Is used as **pain therapy** for treating cramps, pain resulting from injury and chronic pain. The impulses sent to the brain suppress the transmission of pain to the brain via high-frequency impulses and increase of the endorphin release via low-frequency impulses.

Application areas:

- Chronic pain
- Postoperative pain
- Pain due to sports injuries
- Female complaint

EMS: Electrical Muscle Stimulation => Has been used in **sports** and **rehabilitation** medicine for years and it **increases endurance, strengthens muscles and accelerates muscles regeneration**.

With EMS, the electrical impulses stimulate the muscle groups. Muscle fibers contract and relax in a way that mimics natural muscle activation (used for body toning / muscle definition)

Application areas:

- Warming-up of the muscles
- Muscle definition / to tighten muscles and skin
- Muscle relaxation and regeneration (e.g. after marathon)



Electro stimulation (TENS)



MRP: ₹4130/-



MRP: ₹4490/-

EM 27

Ankle

This device can be used for the pain on the ankle and provides pain relief with 4- preprogrammed applications.

- Pain therapy (TENS)
- For use in combatting pain in the ankle
- 2 wear-free contact electrodes made of conductive carbon material
- Ankle cuff (adjustable with hook and loop fastener)
- No contact gel or replacement electrodes required (water contact electrodes)
- 4 pre-programmed applications
- Ankle circumference: approx. 18 - 31 cm

EM 29

Knees & elbows

This TENS device is used for pain therapy for knees and elbows. The snugly, ergonomic shape makes treatment pleasant..

- Pain Therapy
- For use in combating pain in the knee/elbow
- Wear-free water contact electrodes made from conductive carbon material, easy to use
- No contact gel or replacement electrodes required (water contact electrodes)



MRP: ₹5330/-

EM 50

Menstrual Relief

Specially designed for Menstrual pain, this pad has separately adjustable heat and TENS function that can help bring relief and relaxation and has rechargeable batteries with 15 TENS intensity levels.

- Can relieve menstrual or endometriosis pain
- Combination of heat and TENS
- Easy to use and increased flexibility in everyday life with the EM 50
- Reusable
- 2 replaceable gel pads
- 1 heat programme, 15 TENS intensity levels
- LEDs for displaying the set function
- Automatic switch-off after 20 minutes

Electro stimulation

Features	EM 27 Ankle	EM 29 Knees & elbows	EM 50 Menstrual relief
TENS	✓	✓	✓
EMS	X	X	X
Auto switch off	✓ (Countdown timer)	✓ (Countdown timer)	✓ (After 20 min)
Rechargeable batteries	X	X	✓
Intensity level 0- 20	✓	✓	0-15
No of electrodes- 2	✓	✓	✓
4 programs	✓	✓	2



MRP: ₹5245/-



NEW LAUNCH
MRP: ₹11700/-

EM 49

Digital TENS/EMS Unit

This compact device is a combination of various pre-programmed applications namely - TENS, EMS, and massage. These intelligently built simulators can ensure ideal coordination with your therapy routine.

- Pain therapy (TENS), Muscle stimulation (EMS), Relaxation and massage
- Electrode positioning indicator
- Includes 4 self-adhesive gel electrodes (45 x 45 mm)
- 2 separately adjustable channels
- 64 pre-programmed applications
- 6 customisable programs: frequency, pulse width, on/off time programmable
- Adjustable intensity

EM 59

Digital TENS/EMS/HEAT Unit

This digital TENS/EMS is the perfect choice for pain therapy or muscle stimulation. The optional heat function also means it is suitable for relaxation and massage.

- TENS/EMS
- Heat Function
- Electrode positioning Indicator
- Rechargeable battery
- 4 gel pads included
- 3 usage options available:
TENS/EMS+heat, only heat, only TENS
- 2 separately adjustable channels
- 64 pre programmed applications
- 6 customisable programs: frequency, pulse width, on/off time programmable
- Adjustable intensity

FEATURES	EM 49	EM 59
Technology	TENS/EMS/MASSAGE	TENS/EMS/HEAT
Treatment area: Complete Body	✓	✓
No of electrodes: 4	✓	✓
64 Pre-Programmes	✓	✓
Intensity: 0-50	✓	✓
Batteries	3x 1.5V AAA	Rechargeable battery (6x use of 20 min. each)



Air

Modern human beings spend up to 90% of their lives indoors and breathe about 20,000 times (10,000 liters air) a day to supply the body with oxygen. Thus, the indoor climate has a significant impact on health, performance and well-being.

=> Healthy air contributes to:

- Free breathing
- Reduction of dust exposure, especially for allergic persons
- Strengthens the immune system against colds
- Protects skin and eyes from drying out

=> Dry or bad indoor air hinders oxygen transport and oxygen uptake in the respiratory tract and increases the risk for infections and respiratory diseases. Symptoms are e.g., mild discomfort, fatigue and poor concentration but also drying out of furniture, parquet floors, musical instruments, antiques and plants.

=> Dry or bad air can be recognized by e.g., repeated colds and nasal congestion, tightness in the throat, burning eyes, dry skin and chapped lips.





AIR PURIFIERS

An air purifier is a device which removes contaminants from the air in a room to improve indoor air quality. These devices are beneficial for allergy sufferers and asthmatics, and at reducing or eliminating second-hand tobacco smoke..

Functionality:

- The air purifier's fan absorbs the (contaminated) ambient air
- These devices have filters;
 - => The pre-filter filters out larger particles such as dust, dirt and (animal) hairs.
 - => Next, the activated carbon filter filters formaldehyde, benzene and other poisonous gases and odors out of the ambient air.
 - => Lastly, the **EPA / HEPA** filter (depending on the model) filters out small particles such as pollen, bacteria or viruses.
- Additional air cleaning afterwards using either
 - => A selectable ultraviolet light
 - => A selectable ion function
- The cleaned air is emitted via the air outlet back into the room

Result: healthy, clean air!!!



MRP: ₹13,745/-

LR 210

Compact Air Purifier

Clean ambient air while sleeping – with the air purifier with night mode, you can ensure the right indoor environment in your home. Whether it's in your bedroom, living room or other rooms, the filter system gives you clean air by filtering out a large number of airborne viruses and bacteria.

- Air cleaning using a three-layered filter system (pre filter + activated carbon filter + HEPA filter): domestic dust, animal hair, odours, pollen, various bacteria and viruses as well as harmful gases are filtered from the air
- 99.95% of particles $\geq 0.3 \mu\text{m}$ in size are eliminated.
- Additional air cleaning due to the optional ion function
- With night mode: quiet operation with dimmed operating panel
- Illuminated operating panel
- 3 levels
- Timer function
- For rooms up to approx. 26 m^2
- CADR: approx. $100 \text{ m}^3/\text{h}$
- Filter change indicator
- Safety automatic switch-off function



Relaxation and Massage

Massage - what does it do to the body?

When muscle tension, pain or fatigue makes everyday life difficult, a soothing massage can often work wonders. With a massage, the upper skin layers are treated with pressure. The massage promotes blood circulation in the skin and muscles.

In addition, massaging the skin and muscles has a relaxing effect. Hardening in the musculature can be loosened as well. Thus, a massage can also be pain-relieving.

Benefits of Massage ?

Massage improves the general well-being, and in particular:

- Muscular strain
- Tension
- Cellulite
- Blood circulation
- Fatigue symptoms



Relaxation and Massage



MRP: ₹2260/-



MRP: ₹3250/-

MG 21

Infrared Massager

MG 21 massager is designed to provide a soothing and relaxing massage. It comes with three inter-changeable attachments. The device also features an infrared heat mode that can be easily switched-on to intensify the massage whenever needed. It allows you to adjust in between two function levels as per your requirement.

- Soothing and relaxing vibration massage
- Infrared heat can be selected separately
- 3 interchangeable massage attachments (2x nub attachment, brush attachment)
- 2 function levels (vibration massage, vibration massage with infrared heat)

MG 40

Infrared Massager

The device emits a soothing and deep penetrating infrared light which helps in healing. It has 3 function levels namely: vibration massage, infrared heat, and vibration massage with infrared heat. The device helps in relieving muscle tension and pain effectively.

- Soothing vibration massage
- Rotating massage head and 3-way adjustable handle
- Ideal for targeted and spot massage
- Infrared heat can be selected separately
- 4 interchangeable massage attachments (2x nub attachment, brush attachment, cone-shaped attachment)
- 3 function levels (vibration massage, infrared heat, vibration massage with infrared heat)



MRP: ₹6450/-



MRP: ₹7810/-

MG 70

Infrared Massager

MG 70 is designed to offer a powerful tapping massage to keep you relaxed. It comes with two interchangeable attachments—soft attachment and nub attachment. This device also features an infrared heat connectivity, to make your massage session more effective.

- Penetrating tapping massage for relaxation
- Versatile application thanks to removable handle
- Continuously adjustable massage intensity
- Infrared heat can be selected separately.
- 2 interchangeable massage attachments (soft attachment, nub attachment)
- 2 function levels (tapping massage, tapping massage with infrared heat)
- Removable and non-slip handle with soft-touch surface

MG 80


Infrared Massager

This Infrared Massager gives a powerful tap massage to provide ultimate relaxation. It comes up with 2 massage attachments. This massager provides a penetrating massage without applying much force.

- Penetrating tapping massage for relaxation
- Double-head massage
- Continuously adjustable massage intensity
- Infrared heat can be selected separately.
- 2 interchangeable massage attachments (nub attachment, five-point attachment)
- 2 function levels (tapping massage, tapping massage with infrared heat)
- Non-slip handle

Features	MG 21	MG 40	MG 70	MG 80
Vibration Massage	✓	✓	Tapping Massage	Tapping Massage
LED Light/ Infrared heat	✓	✓	✓	✓
Continuously adjustable	X	X	✓	✓
Mains operation	✓	✓	✓	✓
Power	6 W	12 W	22 W	35 W
Attachments	✓	✓	✓	✓





MRP: ₹7800/-

MG 99

Gun Massager

The stylish massage gun can be used to relieve tense muscles in a targeted manner. The massage gun impresses with its lightness and ergonomic shape for comfortable treatment.

- 5 intensity levels
- 4 massage attachments for individual use
- Powerful and soothing trigger point massage
- LED lights for displaying battery level and intensity level.
- Lithium-ion battery for up to 5 hours of use
- Incl. case for practical storage



MRP: ₹3900/-

MG 155

Seat Massager

With its versatile usage and relaxing impact, this vibrating seat cover is definitely an essential for all households! easy to use with its manual control switch.

- For use at home or in the car
- With car adapter for use in the car
- 5 vibration motors
- 2 intensity levels
- 3 individually selectable massage areas
- Heat function with massage or separate
- Easy operation with manual control switch
- Attached using flexible fixing bands
- Can be used on all seats with sufficient seat depth and back rest



MRP: ₹3060/-

FB 12

Footbath Massager

Enjoy a soothing and relaxing foot treatment with the FB12 foot spa. This device performs three main functions. It consumes approximately 60 watts of power. Its portable and light weight nature allows you to carry it anywhere you go. It weights approximately 1.3 kgs.

- With massage attachment
- 3 functions: vibration massage, bubble massage, maintaining water temperature
- Foot rest with massage naps
- Removable splash protection
- Non-slip rubber feet
- Draining outlet



MRP: ₹4285/-

FB 21

Footbath Massager

This foot bath is specifically designed for your exhausted feet after a long tiredness from work. Massager with infrared light field to give a relaxing spa experience.

- With pedicure application
- 3 interchangeable pedicure attachments with handy hook
- 3 functions: vibration massage, bubble massage, maintains water temperature.
- Soothing infrared light field
- Removable roller attachments for foot reflex zone massage
- Foot rest with massage naps
- With practical carrying handle
- Removable splash protection
- Non-slip rubber feet
- Cord winder
- Draining outlet
- Suitable for up to shoe size 49



air
compression
massage



shiatsu
massage



NEW LAUNCH
MRP: ₹17500/-

FM 90

Foot Massager

The Beurer FM90 Shiatsu and Air Compression Foot Massager - the perfect massage for tired feet after a long day. The FM90 works in 3 ways to deliver this standard of massage with a shiatsu massage function, air compression massage function and a heat function.

- Relaxing air pressure massage in 3 intensity levels
- 3 pre-set shiatsu massage treatments
- Optional heat function
- Foot reflexology nodules in the area of the soles and heels
- Removable and washable inner cover
- Suitable up to shoe size 46
- Mains powered: 48 watts.



Light Therapy and Wakeup Lights

What is Light Therapy?

Light is measured in units of "Lux": one lux corresponds to the light of one candle. A human being needs at least 2,500 lux of light per day, but normal electrical light is not sufficient to release sufficient serotonin (a chemical that helps your brain and body communicate and regulate your mood, sleep, appetite, and more).

Even in a well-lit office, e.g., the illumination is just around 500 lux.

Why use a Light Therapy lamp?

- In the medical field, daylight therapy lamps are used as compensation for the effects of a lack of light. Daylight therapy lamps simulate daylight over 10,000 lux. This light can influence the human body and be used as a treatment or as a preventative measure. Normal electric light, however, it is not sufficient to influence the hormonal balance.
- Daylight therapy lamps are intended to compensate the effects of lack of light, particularly sunlight, and provide relief from winter seasonal affective disorders, mood disorders and circadian phase sleep disorders.

The daylight therapy lamps from Beurer simulate daylight over 10,000 lux(distance: ~ 15 - 20 cm)



Light Therapy



MRP: ₹2490/-



MRP: ₹4265/-

IL 11

Infrared Lamp

Beurer IL 11 Infrared Lamp has been developed to treat colds and muscle strains. It comes with an exclusive design and runs on 100 watts. The soothing heat has an analgesic effect on the superficial sensory nerve endings to relieve pain.


- Soothing heat
- For colds and muscle tension
- 5 angle settings
- Max. radiation intensity-680 W/m² at 60 cm

IL 21


Infrared Lamp

Beurer IL 21 Infrared Lamp has been made to treat colds and muscle strains effectively. It comes with a user-friendly design and functions with a 150-watt emitter. The device generates soothing heat that has an analgesic effect on the superficial sensory nerve-endings to relieve pain.

- Soothing heat
- For colds and muscle tension
- 5 angle settings
- Pressed glass bulb
- Medical device



MRP: ₹6280/-



NEW LAUNCH

MRP: ₹9500/-

IL 35

Infrared Lamp

This infrared lamp is made to treat colds and muscle strains. The soothing heat technology of this massager can be used to relieve symptoms and increase well-being. This device has a user-friendly design. It comes with an adjustable lampshade featuring five degrees of inclination.

- Soothing heat
- With 3-level electronic timer
- Switch-off function
- 5 angle settings
- Pressed glass bulb
- Treatment time display in 5-minute intervals
- Max. radiation intensity-550 W/m² at 80 cm

IL 50

Infrared Lamp

The infrared heat lamp is ideal for colds and muscle tension. It has a continuously adjustable shade, a timer with automatic switch off and individually adjustable treatment time.

- Soothing heat
- Ceramic glass plate
- 100% UV Blocker
- Active ventilation
- Overheating protection
- Cord winder

Features	IL 11	IL 21	IL 35	IL 50
Rated power	100 W	150 W	150 W	300 W
Wavelength range/dominant	IRA +IRB / 1100 nm	IRA +IRB / 1100 nm	IRA +IRB / 1100 nm	IRA +IRB / 1250 nm
Protected eyewear included	✓	✓	✓	✓
Angle Settings	5	5	5	Continuously adjustable from 0 – 50°
Switch-off function	X	X	✓	✓
Timer	X	X	3-level electronic timer: 5-10-15 min.	Individually adjustable from 1-15 min.
Replacement bulb (spare part)	✓	✓	✓	✓



Flexible Heating

Providing cosy warmth exactly where it is needed without heating the entire room, therefore, preventing valuable energy from being consumed unnecessarily.

What is the purpose of using an under blanket?

Under blankets keep you warm and comfortable.

It is a perfect solution for back pain or muscle tension. The constant heat radiating from the heated under blankets help reduce pain.

It increases blood flow which in turn releases pressure and any stress accumulated throughout the day.

Beurer Safety System (BSS)

Heated under blankets of Beurer are the newest generation with overheat protection.

- Safety through electronic overheat protection.
- Automatic switch-off in the case of malfunction



Flexible Heating



NEW LAUNCH

MRP: ₹7150/-


HD 75

Cosy Taupe heated over-blanket


This heated over-blanket with fleece fibre surface keeps you nice and cosy on the sofa or in your bedroom.

- Super fleecy surface
- Breathable
- Snugly and skin friendly
- 6 temperature settings
- Beurer Safety System
- Automatic switch off after approx. 3 hrs
- Electronic temperature system
- Illuminated function display
- Removable switch
- Fleece fibres
- Machine washable





MRP: ₹3395/-



MRP: ₹5345/-

HK 42

Super cosy Heat pad

This infrared lamp is made to treat colds and muscle strains. The soothing heat technology of this massager can be used to relieve symptoms and increase well-being. This device has a user-friendly design. It comes with an adjustable lampshade featuring five degrees of inclination.

- Soothing heat
- With 3-level electronic timer
- Switch-off function
- 5 angle settings
- Pressed glass bulb
- Treatment time display in 5-minute intervals
- Max. radiation intensity-550 W/m² at 80 cm

HK 53

Cosy Heat Pad

The provides comforting heat to your neck and back. It is an ergonomically designed electric heating pad. It keeps you warm in winters and helps ease tension particularly in your back or neck. It provides a gentle heat output for longer durations without any inconvenience. The heating pad has an extra soft surface.

- Form fitting design
- Individually adjustable using abdominal belt and soft hook-and-loop fastener at the neck.
- Extra-soft surface
- Breathable
- Fast heating
- 3 temperature settings
- Beurer safety system (BSS)
- Automatic switch-off after approx. 90 min.
- Fleece fibre heating pad
- Illuminated temperature settings.



fitted with
elastic straps



machine
washable



bss safety
system



NEW LAUNCH

MRP: ₹3500/- Single


MRP: ₹5500/- Double



cotton upper
side



machine
washable



3
temperature
settings



MRP: ₹5100/-

TS 15

Under blanket

The under blanket is attached to the mattress, meaning it does not slip and is also able to provide cosy warmth.


- Breathable
- 3 temperature settings
- Beurer Safety System
- Illuminated temperature settings.
- Removable switch
- Top side: fleece; underside: fleece
- Machine washable

UB 33

Under blanket

This under blanket is made with super soft cotton fabric which makes it very comfortable to use.

- Breathable
- 3 temperature settings
- Beurer Safety System
- Illuminated temperature settings.
- Removable switch
- Printed motif
- Top side: cotton; underside: fleece
- Machine washable



MRP: ₹9300/-

UB 56

Under blanket

This breathable mattress is curated with soft and comforting teddy bear fabric. It is XL in size and has two controllers so that the heat can be distinguished on each side. The controller is embedded with LED lights, and there are 4 temperature settings.

- For 2 people
- Breathable
- 2 controllers for 2 individual heat zones
- Top side with soft teddy bear fabric
- 2 zones with 4 temperature settings each
- Automatic switch-off after approx. 12 hours
- Automatic switch-back after 3 hours from level 4/3 to level 2
- Top side: teddy bear fabric; underside: fleece

Features	TS 15	UB 33	UB 56
Breathable	✓	✓	✓
BSS	✓	✓	✓
3 Temperature settings	✓	✓	X
Auto switch off	X	X	✓
Machine washable (30°C)	✓	✓	✓
Material top side- Fleece	✓	Cotton	Teddy bear fabric
Material back side- fleece	✓	✓	✓
Power: 50-60 W	✓	✓	✓
Auto switch back	X	X	Yes, after 3 hours level 4/3 to level 2

Weight and Diagnosis

Weight is an issue of health.

A lot of people are overweight. This is caused by a change in circumstances, undertaking less physical activity, or an excess of high-calorie foods
However, being underweight is also a risk factor

Therefore, regular weight control is important.

Besides weight, a couple of other values can be also determined

1. BMI (Body Mass Index)

It is a measure that uses your height and weight to work out if your weight is healthy. Common method of interpreting the body weight

The weight is set in relation to body height:
$$\text{BMI} = \frac{\text{weight in kg}}{(\text{Height in m}^2)}$$

2. Body fat in %

The human body consists of fat and lean body mass.

A high proportion of body fat is a risk to an increased illness and is not detectable by mere appearance.

3. Body water in %

The human body is made up of **65%** water

The water balance is an important factor for personal well-being: a too low proportion may cause fatigue, cause lack of concentration up to more serious health effects

4. Muscle ratio in %

The musculature is the driving force of your body, well-trained muscles support the bones and relieve the spine.

5. Bone mass

The bone mass of our body is very individual, there is no general guideline

It corresponds to about 12 - 15% of our total weight

Bone mass increases rapidly during childhood and reaches the maximum at about 30 - 40 years in men and women, decreases again with age

Do not confuse bone mass with bone density, which can be determined only at a medical examination (e.g., computed tomography, ultrasound)

Weight and Diagnosis

6. BMR & AMR

BMR (Basal Metabolic Rate) are the calories burned at rest that the body needs to maintain basic functions.

AMR (Active Metabolic Rate) are the calories we consume daily according to our size, age, weight & activity level, while maintaining the current weight.

7. Visceral Fat

There are 2 forms of Fat in our body-

Subcutaneous Fat- Fat that is below the skin and can be felt by touching.

Visceral Fat- Belly fat found deep within the abdominal cavity and surrounds the internal organs, stomach, liver, and intestines. This fat cannot be felt by touching.

Visceral fat is shown in the display as a value between 1 and 59.

Visceral fat percentage value	Rating
1-9	Normal Range
10-14	Elevated Range
15-59	Too high




Weight and Diagnosis



NEW LAUNCH

MRP: ₹3400/-



MRP: ₹3750/-

BF 183

Diagnostic Bathroom Scale

The scale with safety glass and XL Display enables additional measurement of analysis data such as visceral fat and metabolic age.

- Diagnostic bathroom scale with display of visceral fat and metabolic age
- Digit Size: 38 mm
- Standing surface made of safety glass (6mm)
- Weight capacity: 180kg
- Adjustable between: kg/lb/st
- 10 user memories
- Overload indicator, Automatic switch-off
- Modern touch key operation

BF 195

BF 195 Diagnostic Scale has been especially developed for easy weight management. This diagnostic scale has been especially designed to measure the body weight, fat, water, muscle density, bone mass, and calorie requirement (AMR) or the BMI of the user. Additionally, the weighing area comes equipped with a safety glass.

- Weight, body fat, body water, muscle percentage, bone mass, AMR calorie display
- Blue illuminated display
- Modern touch key operation
- 5 activity levels
- User memory spaces: 10
- Switch-on technology: Quick Start
- Digit size: 28 mm
- Automatic switch-off, overload indicator
- Weight capacity: 180 kg

Features	BF 183	BF 195
Displays Visceral Fat & Metabolic Age	✓	X
Illuminated Display	X	✓
Digit Size	38 mm	28 mm
Quick Start	✓	✓
User Memory:10	✓	✓
Batteries	1x3V CR2032	2x3V CR2032







NEW LAUNCH
MRP: ₹4400/-

BF 500 *Bluetooth*

Diagnostic Bathroom Scale

The smart scale enables you to quickly and easily access your measurements any time, anywhere using Beurer Health Manager App.

- Innovative networking between smart phones and scale
- Automatic data transfer within the Bluetooth range after opening the app
- Scale can also be used without the app using the buttons
- With extra large magic display
- Digit height: 40 mm
- Weight, body fat, body water, muscle percentage, bone mass, AMR calorie display
- 8 user memories
- 5 activity level



MRP: ₹2235/-

GS 11

Glass Bathroom Scale

The digital weighing scale has a durable transparent standing surface that is made of high-quality safety glass. It has an elegant design and features smart digital functions. It is extremely easy to use and displays accurate digital measurements.

- Large platform of safety glass
- Easy to read, 25mm digit size
- 150kg capacity
- 100g graduation



MRP: ₹2510/-



MRP: ₹2740/-

GS 10

Personal Bathroom Scale

This bathroom scale has a standing surface that is made of thick safety glass. It features an elegant design with smart digital functions. The device is extremely easy to use and displays accurate digital measurements.

- Easy-to-read LCD display
- Standing surface made of safety glass.
- Flat design height: only 1.9 cm
- Switch-on technology: Quick Start
- Automatic switch-off, overload indicator
- Weight capacity: 180 kg
- Digit size: 26 mm

GS 213

Personal Bathroom Scale

This weighing scale has a wide standing surface that is made of 8mm thick safety glass. The glass is extremely strong and has a high load bearing capacity of up to 180 kgs. It has intelligent automatic functions and elegant design features. It displays extremely accurate measurements on an easy-to-read digital display.

- Large standing surface
- Made of 8 mm safety glass
- Large display
- Flat design height: only 1.8 cm
- Standing surface made of 8 mm thick safety glass
- Switch-on technology: quick start
- Automatic switch-off, overload indicator
- Weight capacity: 180 kg



MRP: ₹3100/-



NEW LAUNCH

MRP: ₹1500/-

PS 160

Personal Bathroom Scale

This weighing Scale has a wide acrylic surface area with an extremely large LCD display. It features smart functions which makes it very efficient in its functionality. It is smartly designed with an accurate weighing mechanism.

- **XXL display**
- Particularly large LCD display
- Standing surface made of comfortable plastic.
- **Decorative element:** glass insert.
- **Switch-on technology:** Quick Start
- Automatic switch-off, overload indicator
- **Weight capacity:** 180 kg
- Digit size: 47 mm

GS 120 *Kompakt Glass*

Personal Bathroom Scale

The Kompakt glass bathroom scale combines space saving design with a modern look and enables precise measurement results.

- **Stylish and compact:** ideal for small spaces
- Elegant mirror effect
- **Digit height:** 22 mm
- **Weight Capacity:** 150 kg
- Overload Indicator



MRP: ₹2950/-

GS 215 *Sans Francisco*

Personal Bathroom Scale

This weighing Scale has a wide acrylic surface area with an extremely large LCD display. It features smart functions which makes it very efficient in its functionality. It is smartly designed with an accurate weighing mechanism.

- **XXL display**
- Particularly large LCD display
- Standing surface made of comfortable plastic.
- **Decorative element:** glass insert.
- **Switch-on technology:** Quick Start
- Automatic switch-off, overload indicator
- **Weight capacity:** 180 kg
- Digit size: 47 mm

Features	GS 10	GS 213	PS 160	GS 120	GS 215
Capacity: 180 kg	✓	✓	✓	150 kg	✓
Graduation: 100 g	✓	✓	✓	✓	✓
Illuminated display	X	X	X	X	✓
Digit Size	26 mm	36 mm	47 mm	22 mm	27 mm
Quick Start	✓	✓	✓	✓	✓
Batteries	1x3V CR2032	2x3V CR2032	1x3V CR2032	1x3V CR2032	3x1.5VCR2032



MRP: ₹2510/-

MS 01

Manual Scale

- Personal scale without battery
- The classic Analogue display of weight
- Available in blue and white
- Easy-to-clean step surface
- Non slippery surface
- 120 kg / 260 lb. capacity
- 100 g / 2 lb. graduation






Babycare



MRP: ₹6145/-



MRP: ₹2510/-

BY 80

Baby Scale

BY 80 is a baby scale that has an automatic and manual hold function which makes it easier to weigh a wriggling baby.


- Automatic and manual hold function for recording the weight of wriggling babies.
- Large LCD display (digit height 23 mm)
- Curved weighing surface
- Switch-on technology: Button-on
- Easy to clean.
- Non-slip, abrasion-resistant rubber feet
- Automatic switch-off

BY 15

Manual Breast Pump

The manual breast pump is ideally suited for occasional pumping or when travelling. It is simple and easy to use. Bottle with anti-colic system for reducing colic and discomfort in infants. With one bottle, manual breast pump attachment, matching screw plug for bottle, bottle stand and a naturally shaped silicone teat for easy feeding incl. cap.

- Bottle with anti-colic system for reducing colic and discomfort in infants.
- Comfortable and gentle due to two pumping levels
- Compact design for easy handling
- Vacuum technology
- Silicone cushion (standard size) with comfortable surface for gentle use
- Accessories can be disinfected



MRP: ₹7980/-

BY 40

Electric Breast Pump

The electric breast pump simulates natural drinking and therefore makes use feel more pleasant. It features a memory function for an optimal personalized mode. Bottle with anti-colic system for reducing colic and discomfort in infants. Very convenient with 10 stimulation levels and 10 pumping level. Easy assembly of parts before use.

- Simulation of natural drinking
- Very convenient with 10 stimulation levels and 10 pumping levels
- Vacuum technology
- Memory function for optimal personalized mode
- Clear and illuminated display.
- Easy assembly of parts before use
- With display of usage time
- Compact design for easy handling
- With a bottle and matching screw plug
- Silicone cushion (standard size) with comfortable surface for gentle use
- Non-slip, abrasion-resistant rubber feet for secure footing
- Accessories can be disinfected.
- Medical device
- Automatic switch-off after 30 min.
- Battery status indicator



Modern And Healthy Lifestyle Solutions For Complete Family By Beurer.



East Zone - Service Center

City	Area	Address	Contact
Assam	Guwahati	INSPIRISYS SOLUTIONS LIMITED H No – 11, 2nd Floor, Triyambak Building, Nabagraha Road, Silpukhuri Guwahati – 781003, Assam, India	Email Id : anwar.hussain@inspirisys.com Phone No : 0361-3518652 Mobile No : 9707639994 Contact Person Name : Anwar Hussain
Bihar	Patna	INSPIRISYS SOLUTIONS LIMITED Office No -110, 1st floor, Surya Prabha Mansions, Off Exhibition Road, Dr. Guha Path Opposite Narmada Apartment, Patna – 800001, Bihar, India	Email Id : sangeeta.srivastava@inspirisys.com Phone No : 0612-3500897,0612-3500897 Mobile No : 9386419001 Contact Person Name : Sangeeta Srivastava
	Muzaffarpur	VISION VIBES TECHNOLOGIES Amber Palace, Jawaharlal Road, Near North Bihar chamber Of Commerce Muzaffarpur – 842001, Bihar, India	Email Id : rajivsingh009@gmail.com Phone No : 7004110809,9334019998 Mobile No : 9334019998 Contact Person Name : Rajiv Kumar
West Bengal	Kolkata	INSPIRISYS SOLUTIONS LIMITED No.19, 1st Floor, R. N. Mukherjee Road, Opposite To Paribahan Bhavan, Kolkata – 700001, West Bengal, India	Email Id : amitdas.s@inspirisys.com Phone No : 033-40033760,033-40033760 Mobile No : 9836646068 Contact Person Name : Amit Das
	Durgapur	EXCEL SERVICES 72 Moulana Azad Sarani City Centre Durgapur – 713216, West Bengal, India	Email Id : excelservicesdgp@gmail.com Phone No : 03432564755 Mobile No : 9933055684 Contact Person Name : Jhumpali Bose
Jharkhand	Jamshedpur	S S ENTERPRISES Shop No.28 2nd Floor Kamani Centre Bistupur Jamshedpur – 831001, Jharkhand, India	Email Id : sse.service@yahoo.com Phone No : 0657-2320821,0657-2320821 Mobile No : 9304564868 Contact Person Name : Sumaghosh
	Ranchi	INSPIRISYS SOLUTIONS LIMITED House No.93, Second Floor, Kadru Maya Bhawan, Old Ag Colony Ranchi – 834002, Jharkhand, India	Email Id : gobinda.chanda@inspirisys.com Phone No : 0651-3501559 Mobile No : 9708503259 Contact Person Name : Gobinda Chanda

	Dhanbad	<p>NAMRATA ENTERPRISES</p> <p>1st Floor, Sahjanand Market Complex, Bank Road</p> <p>Dhanbad – 826001, Jharkhand, India</p>	<p>Email Id : namrata.dhnjhr@gmail.com</p> <p>Phone No : 03262300061,9934990939</p> <p>Mobile No : 9939112404</p> <p>Contact Person Name : Bhupendra Singh</p>
	Ramgarh	<p>SHREYAS EMPIRE</p> <p>NH 33, Bijulia, Beside Reliance Petrol Pump, Opposite KGT Hospital</p> <p>Ramgarh – 829122, Jharkhand, India</p>	<p>Email Id : service.shreyasempire@gmail.com</p> <p>Phone No : 9835516642</p> <p>Mobile No : 7903774006</p> <p>Contact Person Name : Sangeeta</p>
Chhattisgarh	Raipur	<p>RAINA IT HUB</p> <p>Shop No. S-50, Govind Sarang Vyvsaik Parisar, New Rajendra Nagar , Raipur (C.G),</p> <p>Raipur – 492001, Chhattisgarh, India</p>	<p>Email Id : rainaithubryp@gmail.com</p> <p>Phone No : 0771-4098035</p> <p>Mobile No : 7566000011</p> <p>Contact Person Name : Mr. Sanjeev Kumar Chhatri</p>
Odisha	Bhubaneshwar	<p>TECH SOLUTION BHUBANESHWAR</p> <p>Plot No-475, Shaheed Nagar, Infront Of State Bank</p> <p>Bhubaneshwar – 751007, Odisha, India</p>	<p>Email Id : techsolutionbbsr4u@gmail.com</p> <p>Phone No : 0674-2541065,0674-3557621</p> <p>Mobile No : 9437135293</p> <p>Contact Person Name : Guru Krushna Mohanty</p>



West Zone – Service Center

City	Area	Address	Contact
Gujarat	Ahmedabad	INSPIRISYS SOLUTIONS LIMITED 4th Floor, Wing A, B D Patel House, Near Sardar Patel Colony, Naranpura, Ahmedabad – 380014, Gujarat, India	Email Id: kuldeepsinh.rathod@inspirisys.com Phone No : 079-40080022 Mobile No : 9016318499 Contact Person Name: Kuldeep Singh Rathod
	Rajkot	GROWEL INFO TECH Block No. 36, Trisha -2, Astron Society Main Road, Opp. Kings Heights, Amin Marg, Rajkot – 360002, Gujarat, India	Email Id : f1asp.rajkot@gmail.com Phone No : 9824214332,9227776253 Mobile No : 9227776259 Contact Person Name : Satish Kansagara
	Surat	INSPIRISYS SOLUTIONS LIMITED 1st Floor, Flat-H-1+2, Plot -528/A/2, Sagar Plaza, Opp. Yatim Khana Athawa Gate Surat – 395007, Gujarat, India	Email Id : dipika.patel@inspirisys.com Phone No : 6354904288 Mobile No : 9033328120 Contact Person Name : Dipika Patel
	Vadodara	SHREE HARI ELECTRONICS A-107, 1st Floor, Samanvay Silicon, Nr. Dairy Den Circle, Opp Kalyan Hotel, Sayajigunj Vadodara – 390020, Gujarat, India	Email Id : nim_3000@yahoo.com Phone No : 9824331837,9265206575 Mobile No : 9265206575 Contact Person Name : Nimesh Kumar G Patel
		ARIHANT SURGICAL CO Ff-106, Benison Complex, O.P.Road, Near Chali Circle, Opp. Raj Laxmi Society, Vadodara – 390007, Gujarat, India	Email Id : arihantsurgical2002@gmail.com Phone No : 8200990137,7016106315 Mobile No : 9426007427 Contact Person Name : Dipen Patadia
	Mumbai	INSPIRISYS SOLUTIONS LIMITED Shop Nos 16 17 18&19, Vidyarthi Bhuvan, Tribhuvan Road, Near Dreamland Cinema, Grant Road(E), Mumbai – 400004, Maharashtra, India	Email Id : kalpesh.kamble@inspirisys.com Phone No : 022-249764680 Mobile No : 9082778242
Maharashtra		MOBILE CARE SERVICES Office No.40, Ground Floor, Dattani Trade Center, Chandavarkar Road Next To BMC Office, Borivali West Mumbai – 400092, Maharashtra, India	Email Id : mobilcareborivali11@gmail.com Phone No : 9653450133 Mobile No : 7303810756 Contact Person Name : Nisar Ibrahim Siddiqui

	Thane	KRISHNA ELECTRONICS Shop No.1, Sharad Kunj CHS, Backside Of Hdfc Bank, Near Samsung Service Center, Dr. Moose Road, Talav Pali, Thane – 400602, Maharashtra, India	Email Id : krishnathane23@yahoo.com Phone No : 9619562333 Mobile No : 9619562444 Contact Person Name : Anil Kumar Patel
	Yavatmal	STAR ELECTRONICS Shri Balaji Complex, Beside Shyam Talkies, Ram Mandir Road Yavatmal – 445001, Maharashtra, India	Email Id : starelectronicsytl@gmail.com Phone No : 9372721111,8484966264 Mobile No : 8484966264 Contact Person Name : Sanjay G Chhabada
	Pune	INSPIRISYS SOLUTIONS LIMITED 107/2+108, Court Road Office No.202, 2nd Floor, Pasalkar Bhavan, Shivaji Nagar Pune – 411005, Maharashtra, India	Email Id : gajendra.swami@inspirisys.com Phone No : 020-25520041 Mobile No : 9168437937 Contact Person Name : Vishal Pawar Raju Gunte
		NITIN SALES CORPORATION Shop No.1 + Godown, Shree Dattakrupa Society, House No.1411, Sadashiv Peth, Near: Pune Vidhyarthi Gruh Pune – 411030, Maharashtra, India	Email Id : nitinsalespune@gmail.com Phone No : 02024441979,9822097470 Mobile No : 7774932291 Contact Person Name : Akshay Arjun Dhawade
	Navi Mumbai	NANDAI COMPUTERS Office No - 487,C-Wing, Second Floor,Vashi Plaza, Sector-17 Vashi Navi Mumbai – 400703, Maharashtra, India	Email Id : nandaicomputers@gmail.com Phone No : 022-49707338,022-49707338 Mobile No : 8779237660 Contact Person Name : Vishal Shinde
	Nanded	MOBILE DOT STATION Shop No.45, First Floor, Sanman Prestige, Main Railway Station Road, Nanded – 431601, Maharashtra, India	Email Id : mobilestation@gmail.com Phone No : 8432778889 Mobile No : 8983334447 Contact Person Name : Jahangir Mohd. Aslam And Jahangir
	Pimpri	MY MOBILE CARE Umed Bhavan, CTSA No 4715, Canara Bank Building , 2nd Floor, Near Sukhwani, Chamber Station Road , Pimpri Pune Pimpri – 411018, Maharashtra, India	Email Id : mymobilecare2014@gmail.com Phone No : 9890540020 Mobile No : 9890540020 Contact Person Name : Mahesh B Esapurkar

	Nashik	<p>EXCLUSIVE CELLULAR SERVICES</p> <p>Thakkar Bazar</p> <p>Nashik – 422002, Maharashtra, India</p>	<p>Email Id : ecs.nashik@gmail.com</p> <p>Phone No : 7350009893</p> <p>Mobile No : 7350009893</p> <p>Contact Person Name : Sameer Shaikh</p>
	Amravati	<p>AFFLUENT TRADERS</p> <p>2nd Floor Lalita Tower, Near Bheti Hospital, Ambapeth , Rajkamal Chowk,</p> <p>Amravati – 444601, Maharashtra, India</p>	<p>Email Id : afltamt@gmail.com</p> <p>Phone No : 8237681500,8087681500</p> <p>Mobile No : 9890681500</p> <p>Contact Person Name : Ashishb R Bhende</p>
	Kolhapur	<p>ADITYA SOLUTIONS KOLHAPUR</p> <p>No.764, E-Ward, Sudha Shanti Plaza, Gala No.B-5, Near Kanade Hospital, Opposite To Shahupuri 3rd Lane</p> <p>Kolhapur – 416002, Maharashtra, India</p>	<p>Email Id : adityasolutionsklp@gmail.com</p> <p>Phone No : 9763902455,9075082455</p> <p>Mobile No : 9119492455</p> <p>Contact Person Name : Anil Balawant Dhokare</p>
	Andheri West	<p>CELL CARE POINT</p> <p>Shop No 122 1st Floor Pearl Plaza Bldg ,Behind McDonald Opp Andheri Stn</p> <p>Andheri West Mumbai – 400058, Maharashtra, India</p>	<p>Email Id : cellcarepoint98@gmail.com</p> <p>Phone No : 9769992330</p> <p>Mobile No : 9699418605</p> <p>Contact Person Name : Faiz Farooqui</p>
	Nagpur	<p>OM ELECTRONICS</p> <p>Modi Line No.1, Opp. Janki Theatre, Shop No.1, Sitabuldi</p> <p>Nagpur – 440012, Maharashtra, India</p>	<p>Email Id : omelectronicnagpur@gmail.com</p> <p>Phone No : 9028943775</p> <p>Mobile No : 9271232925</p> <p>Contact Person Name : Rupesh Dhomne</p>
	Goregaon East	<p>SPARSH HEALTH CARE</p> <p>D-103 Green Lawn Apartment, Opp: St Plus College, Aarey Road, Goregaon East</p> <p>Mumbai – 400063, Maharashtra, India</p>	<p>Email Id : sparshhealthcare@gmail.com</p> <p>Phone No : 9930086025</p> <p>Mobile No : 9967123446</p> <p>Contact Person Name : Sandeep Shankarrao Suple</p>
	Aurangabad	<p>SHAMS COMPUTERS</p> <p>Shop No.8, Ground Floor, Kalyani's Gurumukh Heights, Cts No.17067,Sheet No.153, H. No.5-1-97, Osmanpura</p> <p>Aurangabad – 431005' Maharashtra' India</p>	<p>Email Id : shamscomputer13@gmail.com</p> <p>Phone No : 9272364898</p> <p>Mobile No : 8552870700</p> <p>Contact Person Name : Khan Iftekhar Ahmed</p>
Goa	Panjim	<p>RAY SERVICES</p> <p>G16-17, B Wing , Jairam Complex, Mala Road, Panjim – 403001, Goa, India</p>	<p>Email Id : rushabhshah@rayservices.co.in</p> <p>Phone No : 7738287583</p> <p>Mobile No : 9604150082</p> <p>Contact Person Name : Rushabh Jayantilal Shah</p>

North Zone - Service Center

City	Area	Address	Contact
Jammu & Kashmir	Srinagar	AL HAYAT HEALTH CARE Hazari Bazar, Near Sheegaz Cinema, Rainawari Srinagar – 190003, Jammu & Kashmir, India	Email Id : alhayathealth@yahoo.com Phone No : 7006189139 Mobile No : 9596250963 Contact Person Name:Wasim Gulzar
Haryana	Karnal	SIMRAN COMMUNICATION Shop No-27, Nirmal Caloni, Behind Mela Ram School Sec-12 Karnal – 132001, Haryana, India	Email Id : rikhi13@gmail.com Phone No : 9215813013 Mobile No : 9416000013 Contact Person Name : Maneet Singh
	Kurukshetra	SHIVOM COMMUNICATIONS Bansal Building, 1st Floor, Near Geeta Dwar, Pipli Chowk Kurukshetra – 136118, Haryana, India	Email Id : shivom.communications5@gmail.com Phone No : 8818045600,8924088880 Mobile No : 8818045600 Contact Person Name : Sanjay & Deepak Bansal
	Gurgaon	SAI OM MOBILE SERVICE Shop No.11, Bus Stand Road, Near-Raghunath Mandir, Opp- Haryana Roadways Work Shop, Gurgaon – 122001, Haryana, India	Email Id : saiom2190@gmail.com Phone No : 9711975393 Mobile No : 9711975393 Contact Person Name : Jitender Kumar
Delhi	Nehru Place	INSPIRISYS SOLUTIONS LIMITED No.601-604, 6th Floor, Ashok Bhavan No 93, Nehru Place New Delhi – 110019, New Delhi, India	Email Id : manas.ghosh@inspirisys.com Phone No : 011-40044880,011-40044880 Mobile No : 1140044880 Contact Person Name : Manas Ghosh
	Wazirpur	AEGIS SOLUTIONS No 331a Rajendera Jaina Tower, Comm. Complex Wazirpur, Wazirpur – 110052, New Delhi, India	Email Id : aegissolutions2016@gmail.com Phone No : 011-47021467 Mobile No : 9999564350 Contact Person Name : Shubham
	Laxmi Nagar	TOUCH AND TYPE A-47 Vijay Block ,Behind Nathu Sweet, Opposite Metro Pillar No.54 Laxmi Nagar New Delhi – 110092,New Delhi, India	Email Id : choudhry.shahid01@gmail.com Phone No : 01140586934 Mobile No : 9354732176 Contact Person Name : Shahid Choudhary

Uttar Pradesh	Lucknow	INSPIRISYS SOLUTIONS LIMITED No.D-78, Nirala Nagar, Near It Crossing Lucknow – 226020, Uttar Pradesh, India	Email Id : madhavi.s@inspirisys.com Phone No : 0522-4953018,8009416234 Mobile No : 8707880021 Contact Person Name : Madhavi Dixit
	Saharanpur	GANPATI INFOTECH 42, Pawan Bela Market, Opp. Income Tax Office, Court Road Saharanpur – 247001,Uttar Pradesh, India	Email Id: ganpatiinfotech@yahoo.co.in Phone No : 9557991334 Mobile No : 9557991334 Contact Person Name : Sanjay Bansal
	Kanpur	MOBILE SERVICE POINT Shop No-5,6 Sagar Planet, Karachi Khana, Kanpur – 208003, Uttar Pradesh, India	Email Id : asheeshkumar784@gmail.com Phone No : 8081373737 Mobile No : 9839088668 Contact Person Name : Asheesh
	Ghaziabad	COMPLETE COMPUTER SOLUTION No.101, 1st Floor, Building No.77,Rahul Complex, Ram Nagar Kirana Mandi Chowk, Ambedkar Road, Ghaziabad – 201001, Uttar Pradesh, India	Email Id : ccsgzb@gmail.com Phone No : 1204313153,9897003153 Mobile No : 9319319100 Contact Person Name : Ritesh K. Agarwal
	Varanasi	NIRMAL INFOTECH 1st Floor, Opp. Union Bank, Shobhnath Dharendra Complex, Rathyatra Varanasi – 221010, Uttar Pradesh, India	Email Id : nirmalinfotech2012@gmail.com Phone No : 9336335060,9336335060 Mobile No : 9336335060 Contact Person Name : Amit
	Gorakhpur	MANSI COMMUNICATIONS Shop No 5, 1st Floor Pratap Market, Cinema Road Golghar Gorakhpur – 273001, Uttar Pradesh, India	Email Id : mansi.caregkp@gmail.com Phone No : 9936370370,9389649444 Mobile No : 9936370370 Contact Person Name : Pradeep Kumar Rai
	Gonda	MOBILE SOLUTION GONDA Opp-Noora Mall Mandir Station Road, Gonda – 271002, Uttar Pradesh, India	Email Id : mobilesolutioun@rediff.com Phone No : 9415534719,7800001001 Mobile No : 9415534719 Contact Person Name : Mr Santosh Thakur
	Jhansi	KRISHNA ENTERPRISES 372/A Janki Complex, Gwalior Road, Civil Line Jhansi – 284001, Uttar Pradesh, India	Email Id : servjhan@gmail.com Phone No : 5102332015,8737922684 Mobile No : 8737922684 Contact Person Name : Sujit Kumar

	Allahabad	<p>CARE AND CARE SERVICES PVT LTD</p> <p>40-C, Lbs Marg,(Opp. A. H. Wheelers Company), Civil Lines</p> <p>Allahabad – 211001, Uttar Pradesh, India</p>	<p>Email Id : anand.careandcare@gmail.com</p> <p>Phone No : 9792970705</p> <p>Mobile No : 9792970705</p> <p>Contact Person Name : Anand</p>
		<p>M R INFOTECH</p> <p>No.5-C, Ambika Tower, Nawab Yusuf Road, Civil Line, Paryagraj</p> <p>Allahabad – 211001, Uttar Pradesh, India</p>	<p>Email Id : manishbharti103@gmail.com</p> <p>Phone No : 9616494185,7860373424</p> <p>Mobile No : 9889769946</p> <p>Contact Person Name : Shivam Kumar</p>
	Shahjahanpur	<p>SHIVA INFOSYS</p> <p>1st Floor, Delhi Tent House, Near Dona & Dona School, Ghanta Ghar</p> <p>Shahjahanpur – 242001, Uttar Pradesh, India</p>	<p>Email Id : shivainfosys1996@gmail.com</p> <p>Phone No : 9235193270,9358554517</p> <p>Mobile No : 9358554517</p> <p>Contact Person Name : Deepak Rathor</p>
	Meerut	<p>S S INFOTECH</p> <p>Sahop No.10 Ashirwas Complex, PI Sharma Meerut</p> <p>Meerut – 250001, Uttar Pradesh, India</p>	<p>Email Id : ssinfotechmrt@gmail.com</p> <p>Phone No : 9997479018</p> <p>Mobile No : 9997479018</p> <p>Contact Person Name : Sachin Kumar Sharma</p>
	Agra	<p>R K INFOTECH</p> <p>Block-33, Shop No. 15b, Chandra Kiran Building, Sanjay Place</p> <p>Agra – 282002, Uttar Pradesh, India</p>	<p>Email Id : rkinfotechagra@gmail.com</p> <p>Phone No : 05624041024,9717766950</p> <p>Mobile No : 9897895024</p> <p>Contact Person Name : Bhanu Kumar</p>
Rajasthan	Udaipur	<p>INSPIRISYS SOLUTIONS LIMITED</p> <p>Office No:329, 3rd Floor University Road, Anand Plaza</p> <p>Udaipur – 313001, Rajasthan, India</p>	<p>Email Id : vinoddutta.sharma@inspirisys.com</p> <p>Phone No : 0294-2429517,9950914279</p> <p>Mobile No : 9950914279</p> <p>Contact Person Name : Vinod Dutt Sharma</p>
	Jaipur	<p>INSPIRISYS SOLUTIONS LIMITED</p> <p>Plot No. A-4, Yudishtar Marg, C-Scheme, Opp Yojana Bhawan,</p> <p>Jaipur – 302001, Rajasthan, India</p>	<p>Email Id : anoop.kumari@inspirisys.com</p> <p>Phone No : 0141-4033605,8009416234</p> <p>Mobile No : 8707880021</p> <p>Contact Person Name : Anoop Kumari</p>
	Alwar	<p>SANDEEP COMPUTER And ELECTRONICS</p> <p>Shah Ji Ka Kuan,Near Jyoti Rao Fule Circle, Bus Stand Road</p> <p>Alwar – 301001, Rajasthan, India</p>	<p>Email Id : sandeepcomputeralwar@gmail.com</p> <p>Phone No : 1442345739,9928016095</p> <p>Mobile No : 9928016094</p> <p>Contact Person Name : NEMI CHAND SAINI</p>

	Jodhpur	VIJAY COMMUNICATIONS Shop No.430, 2nd Floor, Minerva Center,, Mgh Road Jodhpur – 342001, Rajasthan, India	Email Id : vijay.communications2014@gmail.com Phone No : 7568705605,2645825 Mobile No : 9928715825 Contact Person Name : Vijay Prajapat
	Kota	KAY KAY SOFTECH PRIVATE LIMITED No.342-B, Shopping Center, Kota – 324007, Rajasthan, India	Email Id : mahaveer@kaykay.in Phone No : 07442366830,9929909444 Mobile No : 9784590315 Contact Person Name : Mahaveer Sharma
Punjab	Chandigarh	INSPIRISYS SOLUTIONS LIMITED Sco 409 – 410, 2nd Floor, Sector 35 C Chandigarh – 160022, Punjab, India	Email Id : trschd@inspirisys.com Phone No : 0172-5086229 Mobile No : 9501087615 Contact Person Name : Kamal Kishore
	Jalandhar	PERFECT SERVICE POINT 23, Sai Chamber, Ground Floor, Central Market, Jalandhar – 144001, Punjab, India	Email Id : perfectpoint_jal@hotmail.com Phone No : 0181-4631321,0181-4631321 Mobile No : 9417700642 Contact Person Name : Karanjit Singh
	Ludhiana	INNOVATIVE TECHSOFT 3rd Floor, Ranjit Tower, Near Hdfc Bank, Pakhowal, Pakhowal Road Ludhiana – 141002, Punjab, India	Email Id : dlinkcare.ludhiana@gmail.com Phone No : 0161-5006918,0161-5006918 Mobile No : 9781355224 Contact Person Name : Gurupreet Singh
	Patiala	PC GADGETS 59 Factory Area, Near 1st Gol Chakkar (Chanan Chowk) Patiala – 147001, Punjab, India	Email Id : pcgadgets1022019@gmail.com Phone No : 9357880013 Mobile No : 9357880013 Contact Person Name : Jaspreet Singh
	Amritsar	MEDITECH SERVICES 59/1-B, 2nd Floor, Pink Plaza Amritsar – 143001, Punjab, India	Email Id : sharma.anils@yahoo.com Phone No : 9855267715 Mobile No : 9888478535 Contact Person Name : Mr. Anil Sharma
	Pathankot	HINDUSTAN TRADERS S-12/3284 (12), Shiva Ji Market, Opposite Pushap Cinema, Dhangu Road Pathankot – 145001, Punjab, India	Email Id : hindustantraders793@gmail.com Phone No : 7626922588 Mobile No : 9417421958 Contact Person Name : Saksham

South Zone – Service Center

City	Area	Address	Contact
Karnataka	Bangalore	INSPIRISYS SOLUTIONS LIMITED No. 319 To 324, Golden Point 4, 2nd Floor, Off Queens Road, Behind Jasma Devi Bhavan, Vasanth Nagar Bangalore – 560052, Karnataka, India	Email Id : m.gajendran@inspirisys.com Phone No : 080-40870037,080-40870015 Mobile No : 9448689201 Contact Person Name Gajendran
		VSERVE TECHNOLOGIES No 84, 1st Floor, No.1287, 1st Floor, 25th Main Road, 9th Block, Jayanagar, Landmark : Above Poovika Mob Bangalore – 560069, Karnataka, India	Email Id : vservetechnologies.blr@gmail.com Phone No : 080-40997168 Mobile No : 9743232444 Contact Person Name : Mr. Yogesh
	Shimoga	SINGH ELECTRONIC AND COMMUNICATION 3rd Cross, Lbr Complex, 1st Floor Near Shaneshwara Temple, Opp. Baby Care Durgigudi, Shimoga – 577201, Karnataka, India	Email Id : singhcommunication@yahoo.com Phone No : 9844004754,9844066522 Mobile No : 9844004754 Contact Person Name : Satnam Singh
	Belgaum	NEW COSMIC TELECOM Shop:20, Cantonment Complex, Khanapur Road, Camp Belgaum Belgaum – 590001, Karnataka, India	Email Id : newcosmicbelgaum@gmail.com Phone No : 8095525235,8971229580 Mobile No : 9844018837 Contact Person Name : Prakash A Jamaluddin
	Hubli	SHREE ELECTRONICS Behind Pradeep Canteen, Old Income Tax Office Road, Vidyanagar Hubli – 580021, Karnataka, India	Email Id : shree.electronicshubli@gmail.com Phone No : 0836-7968321,9844018837 Mobile No : 9900640561 Contact Person Name : Shivaji
	Mysuru	S K INFOTECH No.108, Ground Floor, KHB Colony, 1st Phase, Kuvempunagar Mysuru – 570023, Karnataka, India	Email Id : skinfotechmys2021@gmail.com Phone No : 0821-4191533,9986389972 Mobile No : 9986538883 Contact Person Name : Kavya P And Shivakumar

	Mangalore	RENUKASHRI ENTERPRISES Vokkaligara Gowdara Bhavana,R. No: 8, 2nd Floor, Chilimbi, Ladyhill Mangaluru – 575006, Karnataka, India	Email Id : sridharamr@outlook.com Phone No : 9448039549 Mobile No : 9449051983 Contact Person Name : Sridhara M R
Andhra Pradesh	Visakhapatnam	INSPIRISYS SOLUTIONS LIMITED Door No.10-50-60/5, Near Raghavendra Swamy Temple, Upstairs Michelin Tyre Shop, TDP Office Road, Waltair Main Road, Ram Nagar Visakhapatnam – 530020, Andhra Pradesh, India	Email Id : madhavi.prasanna@inspirisys.com Phone No : 0891-2513003,0891-2513003 Mobile No : 9533537894 Contact Person Name : Kiran Teja
	Vijayawada	INSPIRISYS SOLUTIONS LIMITED Door No.27-18-44,Opp Maruti Suzuki Spare Parts Congress Office Road , Governor pet Vijayawada – 520002, Andhra Pradesh, India	Email Id : trsvjw@inspirisys.com Phone No : 0866-6663399,0866-6663399 Mobile No : 7671053651 Contact Person Name : Rajeev Ranjan And Sri Manisha
	Rajahmundry	SRI ANNAPURNA MOBILES 8-24-148.Shop.No.5, Vishnu Residence Loge Complex, Main Road Rajahmundry – 533101, Andhra Pradesh, India	Email Id : sriannapurnamobilesrjy@gmail.com Phone No : 0883-6665797,8836665797 Mobile No : 9959686797 Contact Person Name : D.Chinni Krishna
Tamil Nadu	Coimbatore	INSPIRISYS SOLUTIONS LIMITED Designer Complex,, Dr.No -130 B , 2nd Floor, Dr. Nanjappa Road, Coimbatore – 641018, Tamil Nadu, India	Email Id : trscbe@inspirisys.com Phone No : 04224203900 Mobile No : 04224203900 Contact Person Name : Ms. Priyadarshini
	Madurai	INSPIRISYS SOLUTIONS LIMITED 113-B, NCS Complex, Second Floor West Perumal Maistry Street Madurai – 625001, Tamil Nadu, India	Email Id : invmdu_trs@inspirisys.com Phone No : 0452-4208843,0452-4208843 Mobile No : 04524208843 Contact Person Name : Suganya
	Salem	XPLUS COMMUNICATION Shop No:38, Chamundi Complex, Four Roads Salem – 636007, Tamil Nadu, India	Email Id : xplussalem@gmail.com Phone No : 0427-4030771,9080580788 Mobile No : 9080580788 Contact Person Name : Balaji

	Chennai	<p>INSPIRISYS SOLUTIONS LIMITED</p> <p>132,Murugesu Naicker Office Complex Road, IDBI Bank Adjacent Lane, Greams Road, Thousand Light,</p> <p>Chennai – 600006, Tamil Nadu, India</p>	<p>Email Id : saridadevi.r@inspirisys.com</p> <p>Phone No : 044-42015162,044-42015163</p> <p>Mobile No : 9543938810</p> <p>Contact Person Name : Kumaravel.R</p>
	Erode	<p>AARA TECHNOLOGIES</p> <p>No 339/3,4, Arun Complex Brough Road, Near Savitha Hospital Bus Stop</p> <p>Erode – 638001,Tamil Nadu,India</p>	<p>Email Id : sunrisetech.incharge@gmail.com</p> <p>Phone No : 0424-403080,9965568222</p> <p>Mobile No : 9345568222</p> <p>Contact Person Name : Ms.N.Mohanambal And Mr. Joseph</p>
	Hosur	<p>ARS ENTERPRISES</p> <p>No 70k/15a, RCC Complex, Indra Nagar</p> <p>Hosur – 635109, Tamil Nadu, India</p>	<p>Email Id : lenovohosur@gmail.com</p> <p>Phone No : 0434-4244766</p> <p>Mobile No : 9944221332</p> <p>Contact Person Name : Ramar</p>
	Trichy	<p>GAYATHRI COMMUNICATIONS</p> <p>2f-16,Chitra Complex, Near Hotel Raghunath, Chitram Bus Stand</p> <p>Trichy – 620002, Tamil Nadu, India</p>	<p>Email Id : gayathricommunications@gmail.com</p> <p>Phone No : 0431-2708887,0431-2708887</p> <p>Mobile No : 9042698887</p> <p>Contact Person Name : S.Marimuthu</p>
Telangana	Secunderabad	<p>INSPIRISYS SOLUTIONS LIMITED</p> <p>No:32, 1st Floor, Above Canara Bank, Opp: Punjab National Bank, S D Road</p> <p>Secunderabad – 500003, Telangana, India</p>	<p>Email Id : invhyd@inspirisys.com</p> <p>Phone No : 040-27804401</p> <p>Mobile No : 9010454114</p> <p>Contact Person Name : Rahamthunnisa And Srinivas</p>
	Warangal	<p>KAVITHA MOBILES</p> <p>Shop No.08, 1st Floor, Green Square Plaza, Opp: Public Garden Hanumakonda</p> <p>Warangal – 506001, Telangana, India</p>	<p>Email Id : krishnavyaskittu71@gmail.com</p> <p>Phone No : 8978978097</p> <p>Mobile No : 8328339986</p> <p>Contact Person Name : Krishnavyas</p>
Pondicherry	Pondicherry	<p>MAGNANIM SYSTEMS</p> <p>No.2, Bhagath Singh Street, Bharathipuram (Opp To Murugan Steel House)</p> <p>Pondicherry – 605011, Pondicherry, India</p>	<p>Email Id : mail@magnanimsystems.com</p> <p>Phone No : 4134201112,4134201112</p> <p>Mobile No : 9894732828</p> <p>Contact Person Name : Ramesh K</p>

Kerala	Cochin	<p>INSPIRISYS SOLUTIONS LIMITED</p> <p>2nd Floor, Door No.66/1521 (Old No.41/2208), Veekshanam Road, Ernakulam Kochi</p> <p>Cochin – 682018, Kerala, India</p>	<p>Email Id : dileep.purushothaman@inspirisys.com</p> <p>Phone No : 9895702063</p> <p>Mobile No : 9895702063</p> <p>Contact Person Name : Dillep</p>
	Trivandrum	<p>INSPIRISYS SOLUTIONS LIMITED</p> <p>Tc 27/2301(2), First Floor, Lenin Nagar: Rnra-56, Opp.Kmcsu Building Vazhutacadu</p> <p>Trivandrum – 695014, Kerala, India</p>	<p>Email Id : aishwarya.ms@inspirisys.com</p> <p>Phone No : 0471-2324044</p> <p>Mobile No : 9074188523</p> <p>Contact Person Name : Ms. Aishwarya</p>
	Thrissur	<p>REAL VALUE IT SERVICES</p> <p>Tc-13/709/08, 1st Floor, Smart Centre, Near Veterinary Hospital, Kakkala</p> <p>Thrissur – 680021, Kerala, India</p>	<p>Email Id : realvalueitservices@gmail.com</p> <p>Phone No : 8075565583,8075565583</p> <p>Mobile No : 9656141985</p> <p>Contact Person Name : Linto M P</p>
	Kannur	<p>QUBIX TECHNOLOGY</p> <p>Topland Grand Bazar, Bellard Road</p> <p>Kannur – 670001, Kerala, India</p>	<p>Email Id : tech.qubix@gmail.com</p> <p>Phone No : 9497690750,7736301383</p> <p>Mobile No : 9895516922</p> <p>Contact Person Name : Mr. Nandhu And Mr. Shijil A</p>
	Calicut	<p>PRIME INFOTECH</p> <p>No-60/1633, Safa Building, Ram Mohan Road, Chindhavalappu, Near Bairagi Kulam Glp School,</p> <p>Calicut – 673004, Kerala, India</p>	<p>Email Id : sajeevprimeinfo@gmail.com</p> <p>Phone No : 8129823260,9895221251</p> <p>Mobile No : 9656062924</p> <p>Contact Person Name : Jubair</p>

