

VIZBL

Surgical Helmet System
for Joint Replacement Surgeons



VIZBL

Surgical Helmet



Vizbl, the world's lightest surgical helmet system is designed & made-in-India. It is engineered for maximum comfort and AAMI Level 3 protection for the surgeon & patient.



40% Lighter

Vizbl is 40% or 200 gms lighter in comparison to competition. This ensures lower stress on cervical spine during long surgeries.



50% Lower Capital Expenditure

Vizbl eliminates the need for expensive custom Li-ion battery and charging units. Any off-the-shelf 10,000+ mAh power bank, can power Vizbl for 10 hours.



Delivers 540 lpm Airflow

Vizbl delivers 540 liters per min airflow. The airflow in the high performance, Ultra low-noise fans is variable for user comfort.





No Chin-bar design

By eliminating chin-bar, we have maximized surgeon's range of neck movement. The patented visor fixation method creates 190 deg view-angle, enabling peripheral vision.



AAMI Level 3 protection

Vizbl Surgical hood meets AAMI Level 3 protection standards, providing maximum protection to Surgeon & Patient.



Ergonomic design delivers great comfort

Vizbl's ergonomic design ensures great balance & comfort, while retaining the natural Center of Gravity (CoG). Two-way head band controls ensure firm & comfortable seating.



Testimonials

Vizbl Surgical helmet is well balanced and sits well on the head. I was comfortable with the weight and airflow even after 3 hours bilateral surgery.

Dr. Sanjay Londhe, Criticare, Mumbai.

Vizbl Surgical helmet System is immensely quiet and lightweight. These helmets have very good airflow system, keep the surgeon cool and helps them to perform surgeries with ease. I recommend Vizbl to all Orthopedic Surgeons in my network.

Dr. Vividh Makhwana, Apex Hospital, Mumbai.

Vizbl's no-chin bar ergonomic design and significantly lower capital expenses, makes it attractive. Overall its product performance is excellent.

Dr. Sutanu Hazra, AMRI, Mukundapur, Kolkata.



Vizbl is well accepted by Surgeons across India.

9500

Surgeries Performed

180

Helmets in active use



Book your demo, Today!

Scan QR Code to book a demo or email us at pradeep@prodancy.com